

Bard^{*} PerFix^{*} Plug



The Bard^{} PerFix^{*} Plug is designed with pleated edges that conform readily to defects.
The variety of sizes accommodates all groin hernias.*

With more than 4 million implants worldwide, the hernia repair technique using the Bard^{*} PerFix^{*} Plug is proven to be:

Efficient

- Entire operation can take 25 minutes or less²
- Minimal dissection is required
- 4-5cm incision¹
- Local or epidural anesthesia can be used
- Accommodates all groin hernias

Effective

- Utilized in a tension-free repair technique
- Recurrence rates reported at 0.15%¹
- Less than 0.5% chronic pain rate¹

Mesh plug designed for use in a tension-free preperitoneal open inguinal hernia repair technique

The Bard^{*} PerFix^{*} Plug is ideal for use in a tension-free preperitoneal repair technique. Since its introduction in 1993, it has been used in more than four million implants worldwide, and has kept pace with new surgical techniques. For example, in the innovative Modified Technique[†], the outer cone of the Bard^{*} PerFix^{*} Plug, used as an underlay patch, is placed in the preperitoneal plane so that it lies flat against the posterior aspect of the transversalis fascia, assisted by intra-abdominal pressure.

The Bard^{*} PerFix^{*} Plug is designed with pleated edges that conform readily to defects of various sizes and shapes. The inner petals allow the plug to maintain its fluted form and can be removed to customize the Bard^{*} PerFix^{*} Plug to each individual patient. The monofilament polypropylene design ensures healthy tissue ingrowth.

Because surgery with the Bard^{*} PerFix^{*} Plug requires less operating and recovery time compared to conventional hernia repair methods³, procedures can be exceptionally cost-effective.


[†]The choice of technique is subject to the expertise and experience of the trained surgeon.

BARID


DAVOL INC.

TECHNOLOGY
TECHNIQUE
TRAINING
TRUST

MESH PLACEMENT



Classic Technique



Modified Technique

Catalog Number	Quantity	Size	Diameter	
0112750	2/cs.	Small	1.0" x 1.35" (2.5 cm x 3.4 cm)	<input type="checkbox"/>
0112760	2/cs.	Medium	1.3" x 1.55" (3.3 cm x 3.9 cm)	<input type="checkbox"/>
0112770	2/cs.	Large	1.6" x 1.90" (4.1 cm x 4.8 cm)	<input type="checkbox"/>
0112780	2/cs.	Extra Large	1.6" x 2.0" (4.1 cm x 5.0 cm)	<input type="checkbox"/>
0112950	6/cs.	Small	1.0" x 1.35" (2.5 cm x 3.4 cm)	<input type="checkbox"/>
0112960	6/cs.	Medium	1.3" x 1.55" (3.3 cm x 3.9 cm)	<input type="checkbox"/>
0112970	6/cs.	Large	1.6" x 1.90" (4.1 cm x 4.8 cm)	<input type="checkbox"/>
0112980	6/cs.	Extra Large	1.6" x 2.0" (4.1 cm x 5.0 cm)	<input type="checkbox"/>

- I am interested in receiving a procedure introduction kit.
- Please add the Bard® PerFix® Plug to my preference card.
- I would like to have the Bard® PerFix® Plug in stock.
- I am interested in going to a Bard® PerFix® Plug hernia repair training course.

Surgeon's Signature _____

Purchase Order Number _____

Catalog Number _____

Date _____ Quantity _____

References:

- ¹ Millikan KW, Doolas A. A Long-Term Evaluation of the Modified Mesh-Plug Hernioplasty in Over 2,000 Patients. *Hernia*. 2008 June; 12(3): 257-260.
- ² Rutkow, I.M.: The Perfix plug repair for groin hernias. *Surg Clinics of N. America*, 2003; 83; 5: 1079-1098; 78; 6: 1007-1023.
- ³ Kingsnorth Andrew, et al. Lichtenstein Patch or PerFix Plug and Patch in Inguinal Hernia: A Prospective Double-Blind Randomized Controlled Trial of Short-Term Outcome. *Surgery* 2000; 27; 3: 276-283.

*Bard, Davol, PerFix, Ventryo, Visilex, **MK**, Ventralex, Composix, Kugel, Dulex, **CollaMend**, **AlloMax**, **CK**, 3DMax, CruraSoft and PermaSorb are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate. Sepsamesh is a registered trademark of Genzyme Corporation licensed to C. R. Bard, Inc. or an affiliate.

Please consult product labels and inserts for any indications, contraindications, warnings, precautions, and instructions for use.

^{††}Bard® Modified Kugel® Patch

© Copyright 2008, C. R. Bard, Inc. All Rights Reserved.

MMPPSS2

Bard® PerFix® Plug
is just one in a complete family
of hernia repair products:

Ventral Hernia Repair Products

Bard® Ventralex® Hernia Patch
Bard® Composix® L/P Mesh
Ventryo® Hernia Patch
Bard® Sepsamesh® IP Composite
Bard® Dulex® Mesh
Bard® **CollaMend**® Implant
Bard® **AlloMax**® Surgical Graft

Inguinal Hernia Repair Products

Bard® PerFix® Plug
Bard® **MK**^{††} Patch
Bard® **Soft** Mesh
Bard® 3DMax® Mesh
Bard® Visilex® Mesh
Bard® Mesh Flat Sheets and Pre-Shapes

Specialty Products

Bard® **CK**® Parastomal Hernia Patch
Bard® CruraSoft® Patch

Fixation Products

Davol® PermaSorb® Disposable Fixation Device

Bard Surgical Services

Clinical Training Program

National training centers offer instruction in surgical techniques and the ability to view live surgery.

Speaker Program

Educational presentations are given at Grand Rounds, Society Meetings and other venues.

Procedure Introduction Kits

Video programs that describe specific hernia repair techniques and their benefits to you, your patients and your surgical practice.

These services are available for many of the Bard Hernia Repair Products. Please ask your representative, or visit www.davol.com.

BARIED

DAVOL INC.

Davol Inc. • Subsidiary of C. R. Bard, Inc.
100 Crossings Boulevard • Warwick, RI 02886
1.800.556.6275 • www.davol.com
Medical Services & Support 1.800.562.0027